



AB-1650 Super Stat Mini PID



Specifications:

Model Name	AB-1650 Super Stat Mini
Temperature range	-20 to +80°C
Temperature accuracy	±0.05 to ±0.1°C
Temperature display	4 digit LED
Heater	500W SUS spiral sheath heater
Refrigerator	Brushless DC Inverter type 270W Refrigerant R134a
Cooling capacity	550W (473kcal/h : at 24°C) 50Hz
Circulation pump	Maximum flow rate : 4L/min (50Hz), Maximum lift : 3.5m Magnetic sealed pump
Stirring pump	Induction motor / Screw type
Water bath size / Volume	D120 x W150 x H150mm / Approx. 3.5L
Safety function	Leakage overpower breaker, Compressor overheat protection / Circuit Protection, Water level sensor (uppermost and lowermost)
Material	Water bath: SUS, neoprene, Hose: Neoprene, brass Cooling coil: Copper, Heater: SUS Circulation pump: Brass, Float: SUS, Polypropylene
Environment	0 to +30°C, ≤ 85%RH (no condensing)
Size and weight	W192 x D431 x H436mm, Approx.22kg
Standard accessories	Insulated hose (2 pcs), drain hose, user manual



ATTO CORPORATION

3-2-2 Motoasakusa Tel 81 3 5827 4863
Taito-ku Fax 81 3 5827 6647
Tokyo 111-0041 E-mail eig@atto.co.jp
Japan http://www.atto.co.jp

Distributed by: